




THEBES II

Victoreen Model 7020/7040



FEATURES

Flexible and easy-to-use

Range: 50 to 500 cGy/min

No external power is needed at detector location

Dynamic wedge contour mapping

Daily, weekly, and yearly checks as per TG-40

Use with watertank for depth dose calculations

Microsoft® Access database compatible

Microsoft Excel data export

INTRODUCTION


The Victoreen THEBES II consists of a linear ion chamber array, electrometer, communicator, THEBES II Contour Manager software, wall mount power supply, interconnecting cable, and carrying case. An acrylic base plate that holds the ion chamber array and build up plates are also provided.

The linear ion chamber array is permanently connected to the electrometer by a 1.5 m shielded multiconductor cable to eliminate radiation damage to the electrometer.

The THEBES II linear ion chamber array is a waterproof, linear array of 47 ion chambers on 0.5 cm centers. The total active length of the array is 23.42 cm.

The THEBES II communicator provides power and the communication interface for the THEBES II electrometer.

THEBES II Contour Manager Software is a Windows based application that acquires beam contour data from the detector array and display graphically.



60 Degree Dynamic Wedge

APPLICATIONS

The THEBES II is used to perform linear accelerator quality assurance using physicist-preferred air ion chamber technology for dose measurement, avoiding the ion transport issue of the competing liquid chamber technology and radiation damage issues of diode detectors.

The detector array consists of 47 waterproof air ion chambers in a linear array. Two detector arrays are being offered: Model 7020 with 47 ion chambers on a 0.5 cm pitch (23.42 cm total active length) and Model 7040 with 47 ion chambers on a 1 cm pitch (46.88 cm total active length).

The software performs beam contour analysis such as flatness and symmetry. Beam contour data is saved in an Access compatible database, facilitating daily, weekly, and yearly checks, aiding in following the guidelines in TG-40 and TG-51.

Linear array of 47 air ion chambers
0.5 cm spacing
Central axis chamber
Waterproof and vented
Simultaneously checks symmetry,
flatness, light field vs. radiation field
coincidence, field size, beam
center, penumbra, constancy, and
integrated total dose
Real-time dose maps for service
and setup
No electronics in or near the beam
Windows® applications software
Calibrated by Global Calibration
Laboratory (GCL)

SPECIFICATIONS

LINEAR ION CHAMBER ARRAY

Detector Type Ionization chambers, waterproof and vented

Number of Detectors 47

Dimensions

7020 6.2 (w) x 30.2 (d) x 3.7 cm (h)

7040 6.2 (w) x 53.5 (d) x 3.7 cm (h)

Active Area

7020 1.0 x 23.42 cm (23 cm center to center)

7040 0.88 x 46.88 cm (46 cm center to center)

Detector Spacing

7020 0.5 cm

7040 1.0 cm

Inherent Buildup 0.5 cm polystyrene, .02 cm polycarbonate

Inherent Backscatter

Without Mounting Plate 0.3 cm acrylic

With Mounting Plate 0.3 cm acrylic, 1.0 cm acrylic

Radiation Detected

Photons 60 Co to 25 MV

Electrons 6 MeV to 25 MeV

Beam Limits

Maximum Dose/Pulse 12.5 mGy per pulse

Maximum Pulse Rate 1000 pulses per second

Maximum Continuous Dose Rate 500 cGy/min

Ion Chamber Dimensions

7020 0.42 (w) x 0.95 (d) x 0.50 cm (h)

7040 0.88 (w) x 0.88 (d) x 0.50 cm (h)

Ion Chamber Nominal Volume

7020 0.17 cm³

7040 0.36 cm³

Ion Chamber Alignment ± 0.3 mm of chamber outline on top of
array in all axes

Nominal Bias Voltage - 300 V

Weight 0.75 lb (0.34 kg)

ELECTROMETER**Number of Channels** 48**Amplifiers** Non-multiplexed, low leakage, mosfet operational amplifier**Array Scan Time** 1.1 ms**Frame Rate** 5 frames/sec**Status LEDs** 4**Dimensions** 10.9 (w) x 21.6 (d) x 3.2 cm (h)**COMMUNICATOR****Computer Interface** RS-232, DB-9 connector**Power Requirements** 12 VDC, 1 A**Communication Interface****Connector** RJ-45**Baud Rate** 57.6 K**Dimensions** 4.2 (w) x 8.8 (d) x 2.0 cm (h)**SYSTEM****Power Requirements****Input** 120 VAC, 60 Hz, 22 W**Output** 12 VDC, 1 A**Computer Requirements****Computer** IBM® compatible PC, Intel® Pentium® 90 or higher with at least one unused COM port**Operating System** Microsoft Windows 98, ME®, XP®, 2000**Hard Disk Space** 64 MB of available space**Mounting Plate Dimensions****7020** 22.8 (w) x 30.5 (d) x 2.5 cm (t)**7040** 22.8 (w) x 53.8 (d) x 2.5 cm (t)**Calibration****Accuracy** 2%**Reproducibility** 1%**Long Term Stability** 1%**Linearity** 1%**Environmental****Operating Temperature** 50° to 104°F (10° to 40° C)**Storage Temperature** -13° to 149°F (- 25° to 65° C)**Relative Humidity** 20 to 75%, non-condensing**Weight** 21 lb (10 kg)**Accessories Supplied****Wall Mount Power Supply**

Model	Number Description Typical Geo.	Region
14-328	110 VAC 12 VDC 1000 mA	USA, Japan
14-401	230 VAC 12 VDC 1000 mA	Europe
14-414	414 230 VAC 12 VDC 1000 mA	UK
14-414 and 14-416 adaptor	230 VAC 12 VDC 1000 mA	Australia

Model 1090026000 Communication Cable, 75 ft (23 m)

Model 1080024000 Terminator, 120 ohm

Model 1090005000 Communicator

7020

Model 1070006000 Holding Assembly

Model 1070008000 1 cm Buildup Plate

Model 1070007000 2.5 cm Buildup Plate

7040

Model 1070106000 Holding Assembly

Model 1070108000 1 cm Buildup Plate

Model 1070107000 2.5 cm Buildup Plate

Optional Accessories

Model 1070060001 Wellhöfer Watertank Adaptor

7020

Model 1070049000 9.5 cm Buildup Plates

Model 73037020 3-D Watertank Adaptor

Model 73027020 2-D Watertank Adaptor

Model 1070052000 Tray Adaptor Kit, Varian Type 3

Model 1070055000 Tray Adaptor Kit, Elekta

Model 1070061000 Tray Adaptor Kit, Siemens

7040

Model 1070149000 9.5 cm Buildup Plates

Model 73037040 3-D Watertank Adaptor

Model 73027040 2-D Watertank Adaptor

Available Model(s)

Model 7020 THEBES II 20 cm Field Size Array, includes electrometer, communicator, THEBES II Contour Manager software, and carrying case

Model 7040 THEBES II 40 cm Field Size Array, includes electrometer, communicator, THEBES II Contour Manager software, and carrying case

© Elimpex-Medizintechnik, Spechtgasse 32, A-2340 Moedling, Austria
phone +43-2236-410450
fax +43-2236-410459

