



## KVP Meter



The kVp Test-O-Meter family has been specially designed to give accurate readings within seconds.

There are four different models depending on the application:

**RAD/FLU** for measuring all kinds of radiographic and fluoro- scopic equipment.

It automatically calculates kVp down to 10 ms or updates the kVp every second for fluoroscopy.

**MAM** is calibrated for use on a Mo anode with a 30  $\mu\text{m}$  Mo filtration.

**DENT** has a narrow, specially designed detector for measurements on panoramic as well as intra oral equipment.

**CT** with a possibility to adjust for individual calibrations depending on the total filtration of the CT equipment.

---

© Elimpex-Medizintechnik, Spechtgasse 32, A-2340 Moedling, Austria  
phone +43-2236-410450  
fax +43-2236-410459

---

