




## ACR CT Accreditation Phantom Model 468ACR

The new ACR Accreditation phantom model 468ACR contains four modules that test a number of parameters of your CT scanner.

These include:


1. Positioning Accuracy
2. CT # Accuracy.
3. Slice Width
4. Low Contrast Resolution
5. High Contrast (Spatial) Resolution
6. Image Noise

### High Contrast Resolution




- 8 bar resolution patterns.
- Very high contrast compared to background

# High Contrast Resolution




# Low Contrast Resolution

- 1. Contrast is 0.6% from background.
- 2. Groups range from 2mm to 6mm.in Diameter.
- 3. 25mm contrast verification cylinder.



## Linearity and slice width

- 5 Different density areas



Four alignment BBs and Two slice thickness ramps

ACR accreditation gives you the competitive edge. Sign up now at the ACR booth to request your information packet.

---

© Elimpex-Medizintechnik, Spechtgasse 32, A-2340 Moedling, Austria  
phone +43-2236-410450  
fax +43-2236-410459

---

