



Protective Lead Pouch for Storage of LDR Brachytherapy Seeds

Following LDR iodine and palladium seed implantation for prostate cancer, patients may pass seeds during urination. A kit is now available for you to give to your patients to take home for safe storage of the LDR brachytherapy seeds. You no longer have to make the lead pouch, and assemble the components from various locations. This kit consists of a small lead pouch, 0.2mm Pb equivalent, a urologic sieve, a small plastic vial for the seeds, cotton swabs, and patient instructions. The pouch provides safe storage of these LDR brachytherapy seeds and is only 2.75 by 5.5 inches when closed.



Use of the kit helps your medical center comply with 10CFR regulations 35.406 and 20.2003.

Model 90080 contains, 1. P/N 90075 Lead Pouch. 2. P/N 70070 Plastic Vial. 3. P/N70071 Plastic Urological Sieve. Cotton swabs and instructions for the patient are included. When the seeds are in the sieve, they readily adhere to the cotton swabs for transfer to the plastic vial.

Model Number Item

90080	Case of Ten Protective Seed Pouch Kits
	Each Kit Consists of:
	90075 Lead Pouch
	70070 Plastic Vial for Seeds
	70071 Plastic Urological Sieve

This item is available only by the case. There are ten units per case. The minimum order is one case.

© Elimpex-Medizintechnik, Spechtgasse 32, A-2340 Moedling, Austria
phone +43-2236-410450
fax +43-2236-410459

