



MICK® Magazine Sterilizer



The MICK® Magazine Sterilizer, P/N 90087. Magazines are placed in the 12 wells for steam sterilization.

Magazines are stored in 12 wells.

Cover locks in place to secure magazines during transport.

Magazines are completely dried in one drying cycle. A bottom vent allows water to drain.

Thick aluminum walls around the magazines provides maximum shielding.

More convenient than the lead shipping containers.

MICK® is a registered trademark of Mick Radio-Nuclear Instruments, Inc.

Doc. No. 1145-00, One page, 08.03.99

© Elimpex-Medizintechnik, Spechtgasse 32, A-2340 Moedling, Austria
phone +43-2236-410450
fax +43-2236-410459

