



Seed Sterilization Pill Box

For sterilization of Iodine and Palladium seeds

Cover can be secured for safe handling
Rounded inner corners keep seeds positioned for easy removal
Top flips over for secondary receptacle
Solid construction assures stability



The Seed Sterilization Pill Box, Part Number 90086, is designed to allow sterilization of a large number of iodine or palladium seeds at one time. The seeds are placed in the bottom of the well and covered with the screen. The cover can be secured for safe handling via the interlocking pins. With the cover locked, the Seed Sterilization Pill Box can be safely carried to the autoclave. To sterilize, remove the cover, place the Seed Sterilization Pill Box with the screen in place in the autoclave and sterilize. Stainless steel construction, fully autoclavable. Steam autoclave only. Height: 1.5 in, 3.8 cm; Diameter: 2.5 in, 6.4 cm; Walls: 3/8"; Weight: 1.6 lb., 0.73 kg. Designed in collaboration with J. Larry Beach, PhD, Lexington Radiation Therapy Center, Lexington, KY

© Elimpex-Medizintechnik, Spechtgasse 32, A-2340 Moedling, Austria
phone +43-2236-410450
fax +43-2236-410459

