

# Delta4 Family

Choose your patient QA tools to fit your process

## Delta4 Phantom+ MR

### PRE-TREATMENT

#### Accurate patient QA in less time than ever

The Delta4 Phantom+ MR is a specialized version of the Delta4 Phantom+ designed for use with MR-linacs. Sharing the same qualities as the original phantom, the Delta Phantom+ MR is MR-compliant for field strengths up to 1.5T and offers the most efficient and accurate verification of Elekta Unity and ViewRay MRIdian.



it an ideal tool for characterizing complex treatment fields

- Compatible with multiple treatment machines, and can be integrated seamlessly into existing patient QA workflows, ensuring minimal disruption to your clinic's operations, and providing confidence in the delivery of accurate treatments.



- The ideal solution for patient QA in radiation therapy for MRI-linacs
- Unique design allows for accurate and reliable measurements of dose distributions in MR-linac environment
- Equipped with high-resolution detector array that enables sub-millimeter spatial resolution, making